

# CLEANING & CARE

## DENTAL HANDPIECES

veterinary technology for life  
**Eickemeyer®**

*For maximum durability.*

*With regular care and cleaning after each procedure, you can maintain the full functionality of your handpieces.*

1) After each procedure, remove the handpiece from the coupling and then remove the drill (ideally with a pair of tweezers).

2) **External cleaning of the handpieces/turbines**

The handpiece should be cleaned with warm, running water (e.g. with a toothbrush) so that abrasive dust or saliva is washed away. Never clean the handpieces in an ultrasonic bath!

3) **Disinfection of handpieces/turbines**

The handpieces can be treated externally in the spray water/spray air channels with disinfectants. Never immerse handpieces in disinfectant solution!

4) **Care of handpieces/turbines**

Handpieces should be sprayed with care oil after each clean. Spray the care oil into the water and spray channels of the handpieces for 1 – 2 seconds. If the liquid emerging at the head is still mixed with dirt, the process should be repeated.



*Injecting maintenance oil into the spray channel (maintenance oil item no. 174010 and spray channel adapter item. no. 17401150)*



*Injecting care oil into handpieces with a planar receiving shaft (care oil item no. 174010 and handpiece adapter item no. 17401155)*

5) **Care of collet chucks**

The free movement of the push buttons is very important for the fixed clamping function of the drill. The clamping forceps must be sprayed with the maintenance oil at least once a week from the drill side. Only use diamond grinders and drills from well-known manufacturers. If the shaft material is too soft, scrapings may occur if the drill is twisted in the clamping jaw. This results in fine chips, which can become caught in the clamping system. Thereafter, the drill may no longer be removed and the tensioning system thereby damaged.

6) **Care**

It is important to keep the all instrument joints well lubricated with instrument oil (e.g. no. 563707) to avoid any friction or corrosion. Other care checks include smoothness of movement. We offer a repair service for instruments.

“  
**OIL REGULARLY FOR  
A LONG LIFESPAN**

”

**EICKEMEYER®**

**GERMANY** | Eltstraße 8 | 78532 Tuttlingen | Deutschland | T +49 7461 96 580 0 | F +49 7461 96 580 90 | info@eickemeyer.de | www.eickemeyer.de

**CANADA** | 250 Briarhill Dr. | Stratford | Ont. Canada | N5A 7S2 | T +1 519 273 5558 | F +1 780 963 5428 | info@eickemeyervet.com | www.eickemeyervet.com

**DENMARK** | Lysbjergvej 6, Hammelev | 6500 Vojens | Denmark | T +45 7020 5019 | F +45 7353 5019 | info@eickemeyer.dk | www.eickemeyer.dk

**ITALY** | Via G. Verdi, 8 | 65015 Montesilvano (PE) | Italia | T +39 085 935 4078 | F +39 085 935 9471 | info@eickemeyer.it | www.eickemeyer.it

**NETHERLANDS** | Bedrijventerrein Pavijen-West | Bellweg 44 | 4104 BJ Culemborg | Nederland | T +31 345 58 9400 | info@eickemeyer.nl | www.eickemeyer.nl

**POLAND** | AL Jana Pawła II 27 | 00-867 Warszawa | Polska | T +48 22 185 55 76 | F +48 22 185 59 40 | info@eickemeyer.pl | www.eickemeyer.pl

**SWITZERLAND** | Sandgrube 29 | 9050 Appenzell | Schweiz | T +41 71 788 23 13 | F +41 71 788 23 14 | info@eickemeyer.ch | www.eickemeyer.ch

**UNITED KINGDOM** | 3 Windmill Business Village | Brooklands Close | Sunbury-on-Thames | Surrey | TW16 7DY | UK | T +44 20 8891 2007 | www.eickemeyer.co.uk