



EICKEMEYER® Aspheric  
Ophthalmoscopy

veterinary technology for life  
**Eickemeyer**®

# *Look into my eyes...*



173474



173475

## **EICKEMEYER® Aspheric Ophthalmoscopy**

The 20 D and 30 D lenses have been known as the industry standard for decades. They are the most commonly used for general diagnosis/treatment.

The large diopter gives a small image for a good overview. Small diopter gives you a bigger image with less overview. The bigger the eye, the „smaller“ the diopter of the lens needed.

► Cat/Dog: 20 D or 30 D

20 D  
**173474**

30 D  
**173475**

**EICKEMEYER®**  
Eltastraße 8  
78532 Tuttlingen  
Germany

T +49 7461 96580 0  
F +49 7461 96580 90  
export@eickemeyer.com

**TELEPHONE +49 7461 96580 0**  
**www.eickemeyer.com**