



EickWrap FUN
Cohesive Bandages

veterinary technology for life
Eickemeyer®

Colour up your life



EickWrap FUN Cohesive Bandages

Introducing EickWrap Bandages with fun, colourful motifs!
The tried and tested EickWrap bandages are now available in vibrant designs. EickWrap's self-adhesive surface ensures excellent cohesion without sticking to skin, hair, or other textiles. These bandages do not constrict blood circulation or joint movements, making them ideal for dog and cat paws, as well as horse legs and hooves.

- ▶ 6 Colourful Motifs: Choose from a variety of eye-catching designs
- ▶ Highest Quality Cohesive Bandage: Reliable and durable
- ▶ Conforms to Shape: Adapts to the area where it is applied
- ▶ Breathable: Allows air circulation over the skin and around the dressing
- ▶ Individually Packaged: Convenient and hygienic
- ▶ Contains Latex
- ▶ Length: 4.5 metres per roll
- ▶ Width: 5 cm
- ▶ Pack of 36 Pieces

050758

EICKEMEYER®
Eltastraße 8
78532 Tuttlingen
Germany

T +49 7461 96 580 0
F +49 7461 96 580 90
export@eickemeyer.com

TELEPHONE +49 7461 96 580 0
www.eickemeyer.com