



EICKEMEYER® LED  
Wood's Lamp

veterinary technology for life  
**Eickemeyer®**

# *Take a closer look*



## ***EICKEMEYER® LED Wood's Lamp***

The Wood's Lamp is a special diagnostic tool designed to illuminate and magnify skin, enabling clinicians to accurately inspect and diagnose various conditions. This lightweight and user-friendly LED lamp helps identify fluorescent diseased areas including bacterial infections, and pigment disorders.

- ▶ Window diameter: 65 mm
- ▶ 1.5x magnification
- ▶ Working distance: approx. 10 cm
- ▶ High depth of field
- ▶ 16 LEDs for uniform illumination
- ▶ Wavelength: 365 nm
- ▶ Weight: 305 g

**179030**

**EICKEMEYER®**  
Eltastraße 8  
78532 Tuttlingen  
Germany

T +49 7461 96 580 0  
F +49 7461 96 580 90  
[export@eickemeyer.com](mailto:export@eickemeyer.com)

**TELEPHONE +49 7461 96 580 0**  
[www.eickemeyer.com](http://www.eickemeyer.com)