

veterinary technology for life

# Eickemeyer®

[www.eickemeyer.com](http://www.eickemeyer.com)

Telephone +49 7461 96 580 0

## *EickLoxx Small Contour Plate Set*




# EICKLOXX SMALL CONTOUR PLATE SET – COMPONENTS AND CHARACTERISTICS

The EickLoxx Small Contour plate set includes various anatomically preformed osteosynthesis plates for the treatment of fractures of the distal radius and femur, ilium, acetabulum and for arthrodesis of the carpus and tarsus. It complements the cut-to-size plates of the EickLoxx Small osteosynthesis system and is therefore also designed for the treatment of cats and dogs up to 15 kg. The instruments from the EickLoxx Small osteosynthesis system are used in a similar way.

The titanium plates of the EickLoxx Small Contour plate set are available in two plate thicknesses: The 2.2 mm plates (gold, green and magenta) can be polyaxially locked with the 1.7 mm and 2.3 mm screws, and the 3.0 mm plates (light blue) can be polyaxially locked with the 2.7 mm, 3.5 mm and 4.0 mm screws from the EickLoxx series.

Despite contouring, the plates can be bent as required using the bending irons and / or the plate bending pliers with rollers and shortened using the cutting pliers.

The plate tray with the EickLoxx Small Contour plates fits into the container (Item No. 185555) on the EickLoxx Small Instrument and Implant Tray.



## Titanium EickLoxx Small Contour Bone Plates

- ▶ Titanium alloy for maximum biocompatibility
- ▶ Multidirectional locking bi- and monocortical
- ▶ Eliminates abrasion
- ▶ Geometry optimised strength
- ▶ Bendable in two planes



**197910 – 197915, 197920 – 197925, 197930 – 197935**

## Characteristics

- ▶ EickLoxx bone plates minimise contact with the periosteum and thus reduce the iatrogenic stress on bone perfusion that is common with conventional compression plates
- ▶ Preservation of perfusion significantly reduces the risk of infection and accelerates bone healing
- ▶ Resistance to infection is also increased by the biocompatibility of titanium and the absence of friction wear

## Biological advantages

- ▶ Reduces damage to the vascularisation
- ▶ Increased resistance to infection
- ▶ Accelerated healing



## EICKLOXX SMALL CONTOUR PLATE SET – ITEM LIST

EickLoxx Small Contour Plate Set		
Item no.	Description	Quantity
197900	Complete set, consisting of:	
197950	EickLoxx Small Contour Plate Tray, without implants	1
197910	EickLoxx Small Contour Bone Plate, L-plate left, 9 holes, system 1.7 / 2.3 mm, gold, dimensions (in mm): L 39 x W 18 x H 2.2	1
197911	EickLoxx Small Contour Bone Plate, L-plate right, 9 holes, system 1.7 / 2.3 mm, gold, dimensions (in mm): L 39 x W 18 x H 2.2	1
197912	EickLoxx Small Contour Bone Plate, T-plate, 10 holes system 1.7 / 2.3 mm, gold, dimensions (in mm): L 54 x W 18 x H 2.2	1
197913	EickLoxx Small Contour Bone Plate, Y-plate, 10 holes, system 1.7 / 2.3 mm, gold, dimensions (in mm): L 56 x W 18 x H 2.2	1
197914	EickLoxx Small Contour Bone Plate, Hockey-plate left, 9 holes, system 1.7 / 2.3 mm, gold, dimensions (in mm): L 58 x W 11 x H 2.2	1
197915	EickLoxx Small Contour Bone Plate, Hockey-plate right, 9 holes, system 1.7 / 2.3 mm, gold, dimensions (in mm): L 58 x W 11 x H 2.2	1
197920	EickLoxx Small Contour Bone Plate, C-plate, 4 holes, system 1.7 / 2.3 mm, green, dimensions (in mm): L 28 x W 6 x H 2.2	1
197921	EickLoxx Small Contour Bone Plate, C-plate, 5 holes, system 1.7 / 2.3 mm, green, dimensions (in mm): L 30 x W 6 x H 2.2	1
197922	EickLoxx Small Contour Bone Plate, C-plate, 6 holes, system 1.7 / 2.3 mm, green, dimensions (in mm): L 40 x W 6 x H 2.2	1
197923	EickLoxx Small Contour Bone Plate, Cricket-plate, 4+4 holes, system 1.7 / 2.3 mm, magenta, dimensions (in mm): L 55 x W 6–8 x H 2.2	1
197924	EickLoxx Small Contour Bone Plate, Cricket-plate, 6+6 holes, system 1.7 / 2.3 mm, magenta, dimensions (in mm): L 80 x W 6–8 x H 2.2	1
197925	EickLoxx Small Contour Bone Plate, Cricket-plate, 8+8 holes, system 1.7 / 2.3 mm, magenta, dimensions (in mm): L 105 x W 6–8 x H 2.2	1
197930	EickLoxx Small Contour Bone Plate, L-plate left, 9 holes, system 2.7 / 3.5 mm, light blue, dimensions (in mm): L 56 x W 24 x H 3.0	1
197931	EickLoxx Small Contour Bone Plate, L-plate right, 9 holes, system 2.7 / 3.5 mm, light blue, dimensions (in mm): L 56 x W 24 x H 3.0	1
197932	EickLoxx Small Contour Bone Plate, T-plate, 9 holes, system 2.7 / 3.5 mm, light blue, dimensions (in mm): L 72 x W 22 x H 3.0	1
197933	EickLoxx Small Contour Bone Plate, Y-plate, 9 holes, system 2.7 / 3.5 mm, light blue, dimensions (in mm): L 72 x W 18 x H 3.0	1
197934	EickLoxx Small Contour Bone Plate, Hockey-plate left, 9 holes, system 2.7 / 3.5 mm, light blue, dimensions (in mm): L 75 x W 15 x H 3.0	1
197935	EickLoxx Small Contour Bone Plate, Hockey-plate right, 9 holes, system 2.7 / 3.5 mm, light blue, dimensions (in mm): L 75 x W 15 x H 3.0	1

EICKEMEYER KG  
Eltastraße 8  
78532 Tuttlingen  
Germany  
T +49 7461 96 580 0  
F +49 7461 96 580 90  
export@eickemeyer.com  
www.eickemeyer.com