



X-Ray
Protective Goggles

veterinary technology for life
Eickemeyer®

For your eyes only



X-Ray Protective Goggles

The lens of the eye is particularly sensitive to radiation. Therefore, it is essential for operators of X-Ray equipment to wear protective glasses to shield their eyes from harmful scatter radiation and prevent long-term eye damage. A well-fitting pair of protective glasses provides not only frontal but also lateral protection.

- ▶ High-quality plastic frame
- ▶ Pb front protection glass: 0.75 mm
- ▶ Pb side protection: 0.5 mm
- ▶ Goggle strap with width adjustment
- ▶ Incl. case

701010