

Veterinary Consumables Catalogue

pro
veterinary experts



EICKINJECT SYRINGES

- ▶ 3 part syringe with Luer slip connector
- ▶ Individually packaged – sterile

- ▶ 1 ml, with displacement spike, centric
- ▶ Pack of 100 pieces
- 050356**

- ▶ 2 ml, centric
- ▶ Pack of 100 pieces
- 050357**

- ▶ 3 ml, centric
- ▶ Pack of 100 pieces
- 050358**

- ▶ 5 ml, centric
- ▶ Pack of 100 pieces
- 050359**

- ▶ 10 ml, eccentric
- ▶ Pack of 100 pieces
- 050360**

- ▶ 20 ml, eccentric
- ▶ Pack of 50 pieces
- 050361**

- ▶ 50 ml, eccentric
- ▶ Pack of 25 pieces
- 050362**



050356

050357

050358

EICKINJECT HYPODERMIC NEEDLES

- ▶ Individually packaged – sterile
- ▶ Pack of 100 pieces

- ▶ Orange
- ▶ 25G x 5/8"
- 050376**

- ▶ Orange
- ▶ 25G x 1"
- 050369**

- ▶ Blue
- ▶ 23G x 5/8"
- 050372**

- ▶ Blue
- ▶ 23G x 1"
- 050368**

- ▶ Black
- ▶ 22G x 1 1/4"
- 050367**

- ▶ Green
- ▶ 21G x 5/8"
- 050371**

- ▶ Green
- ▶ 21G x 1"
- 050370**

- ▶ Green
- ▶ 21G x 1 1/2"
- 050366**

- ▶ Yellow
- ▶ 20G x 1"
- 050375**

- ▶ Yellow
- ▶ 20G x 1 1/2"
- 050365**

- ▶ Ivory
- ▶ 19G x 1 1/2"
- 050364**

- ▶ Pink
- ▶ 18G x 1 1/2"
- 050363**



050376

050369

050372

050368

050367

050371



050370

050366

050375

050365

050364

050363

I.V. GIVING SET

- ▶ Pressure / gravity infusion
- ▶ Individually packaged sterile sets
- ▶ DEHP & Latex free
- ▶ Compatible with infusion pump InfusoVet, InfusoVet Easy II, InfusoVet Connect and InfusoVet IP
- ▶ Pack of 100 pieces
- ▶ Extension 050166 is compatible with this giving set
- ▶ Please note: If you use a pump that is not from EICKEMEYER®, your pump for this infusion set must be calibrated accordingly

402071

**ORDER
ONLINE**

www.eickemeyer.com



I.V. GIVING SET FLOW REGULATOR

- ▶ Flow rate: 5 – 250 ml per hour

051087



I.V. GIVING SET

- ▶ 20 drops/ml
- ▶ Length: 1.50 m
- ▶ Y-connector and Filter
- ▶ Includes Luer-Lock connector
- ▶ Individually packaged – sterile
- ▶ Pack of 50 pieces

050165



INFUSION SET FOR HORSES

- ▶ Easy removal by large volume hoses
- ▶ Max. Flow rate (with drip control): 0.15 l / min
- ▶ Max. Flow rate (without drip control): 0.5 l / min
- ▶ Tubing length: 505 cm
- ▶ DEHP-free
- ▶ Latex-free
- ▶ Sterile

402135



I.V. GIVING SET EXTENSION

- ▶ Includes Luer-Lock connectors
- ▶ Length: 1.50 m
- ▶ Also compatible with giving set 402071
- ▶ Individually packaged – sterile
- ▶ Pack of 50 pieces

050166



I.V. GIVING SET EXTENSION

- ▶ Straight
- ▶ Length: 140 cm

050161



SURGICAL DRAPES, PVC

PVC

- ▶ Made of soft, tear-resistant PVC
- ▶ Excellent adhesive properties with spray adhesive
- ▶ Individually packaged – sterile

- ▶ Dimensions (in cm):
L 30 x W 45
- ▶ Pack of 25 pieces

920145

- ▶ Dimensions (in cm):
L 120 x W 120
- ▶ Pack of 10 pieces

921212

- ▶ Dimensions (in cm):
L 60 x W 90
- ▶ Pack of 25 pieces

920190

- ▶ Dimensions (in cm):
L 120 x W 180
- ▶ Pack of 10 pieces

921218

- ▶ Dimensions (in cm):
L 90 x W 120
- ▶ Pack of 10 pieces

920212

- ▶ Dimensions (in cm):
L 120 x W 250
- ▶ Pack of 10 pieces

921225

Coated PVC

- ▶ PVC coating with absorbent layer
- ▶ Excellent adhesive properties with spray adhesive
- ▶ Individually packaged – sterile

- ▶ Dimensions (in cm):
L 35 x B 50

921500

- ▶ Dimensions (in cm):
L 50 x W 50

921501

- ▶ Dimensions (in cm):
L 50 x W 75

921502

- ▶ Dimensions (in cm):
L 75 x W 90

921503

- ▶ Dimensions (in cm):
L 100 x W 150

921504

- ▶ Dimensions (in cm):
L 150 x W 200

921505



920190

3M™ IOBAN™ ANTIMICROBIAL INCISE DRAPES

The Ioban™ antimicrobial incision drapes are designed to reduce the risk of post-operative wound infections. They are applied to the patient's skin at the incision site in order to create a sterile operating field and to develop an antimicrobial effect during the surgical procedure. Ioban™ is an iodophor-impregnated, adhesive incision film that forms a barrier and has an antimicrobial spectrum of activity on the patient's skin flora.

952600

Size: Small

- ▶ Dimensions (in mm):
L 150 x W 200
- ▶ Adhesive Area (in mm):
L 100 x W 200

191600

Size: Medium

- ▶ Dimensions (in mm):
L 600 x W 350
- ▶ Adhesive Area (in mm):
L 350 x W 350

191601



EICKSKIN ADHESIVE SPRAY

- ▶ Quick drying adhesive, which is easy and convenient to apply
- ▶ Suitable to adhere drape(s) to skin
- ▶ Content: 200 ml

952600



NASOGASTRIC FEEDING TUBE FOR DOGS

**ORDER
ONLINE**

www.eickemeyer.com

- ▶ Suitable for canines
- ▶ PVC construction
- ▶ Closed, rounded tip with 2 eyelets
- ▶ Can also be used as a urinary catheter for cattle

- | | |
|------------------|------------------|
| ▶ Ch. 18 | ▶ Ch. 25 |
| ▶ Ø 6 mm | ▶ Ø 8 mm |
| ▶ Length: 800 mm | ▶ Length: 800 mm |

460006

460008



460006



460008

GASTRIC LAVAGE SET INC. 3 PROBES

This gastric lavage set includes:

- ▶ Fluid container
- ▶ 3 probes : 10 mm, 20 mm, 30 mm diameter

460620



REDON DRAIN 6 FG

- ▶ Even circular perforation provides numerous outlets for drainage
- ▶ Radiopaque stripes
- ▶ Individually packaged - sterile

456106



MIN-REDON BELLOWS

- ▶ Low vacuum system with active suction for use in small animals
- ▶ Universal screw connector for 6 Fg and 8 Fg Redon drains
- ▶ Volume: 20 ml
- ▶ Individually packaged - sterile

457570



MINI REDOVAC DRAIN SET 6 FG

- ▶ Low vacuum wound drainage system with active suction
- ▶ Includes a drain with an attached needle
- ▶ Bellow Volume: 20 ml
- ▶ Tube Length: 35 cm

456006



CHEST TUBE VALVE ONE-WAY

- ▶ Ideal device for restoring the negative pressure in the pleural cavity and for general suction drainage
- ▶ PVC construction
- ▶ LEO valve

457580



PENROSE DRAIN

- ▶ Made of silicone
- ▶ Length: 100 cm
- ▶ Individually packed - sterile
- ▶ Pack of 10 pieces



Size 1

- ▶ Ø 6.3 mm
- ▶ Width: 10 mm (flat)

457001

Size 3

- ▶ Ø 10.0 mm
- ▶ Width: 16 mm (flat)

457003

Size 2


- ▶ Ø 8.0 mm
- ▶ Width: 13 mm (flat)

457002

Size 4

- ▶ Ø 13.0 mm
- ▶ Width: 20 mm (flat)

457004



URINARY CATHETER CANINES (MALE)

- ▶ Suitable for canines (male)
- ▶ Rounded tip
- ▶ Two lateral eyes
- ▶ Plastic construction
- ▶ Luer cone with closing cap
- ▶ Packaged/sterile
- ▶ L 500 mm

- ▶ Grey
- ▶ Ø: 1.5 mm

451215

- ▶ Red
- ▶ Ø: 2.5 mm


451225

- ▶ Green
- ▶ Ø: 2.1 mm

451221

- ▶ Black (CH 10)
- ▶ Latex free
- ▶ Ø: 3.3 mm

451233



KATKATH® CATHETER 3.5FR

- ▶ Allows passage of urine when the urethra has a functional or physical obstruction
- ▶ Constructed from non-reactive material which has similar properties to Teflon but is more flexible
- ▶ Kink and crush-resistant so no stylet required
- ▶ Suitable for cats

Features:

- ▶ Non-reactive material
- ▶ Low coefficient of friction
- ▶ Lightweight
- ▶ Radiopaque
- ▶ Kink & Crush-resistant

Size 1

- ▶ L 110 mm

451100

Size 3

- ▶ Ø 1.1 mm

- ▶ L 180 mm

- ▶ Pack of 10 pieces

451103

Size 2

- ▶ L 140 mm

451101



URINARY CATHETER FOR TOMCAT

Before use, add 3 ml of saline or sterile water to the package to activate the hydrophilic coating.

- ▶ Cylindrical, rounded tip
- ▶ 2 lateral eyes
- ▶ With Luer-Lock cone
- ▶ With stylet

- ▶ Individually sterile packaged
- ▶ Pack of 12
- ▶ 3 Fr x 5 1/5"

- ▶ Length: 130 mm
- ▶ Ø: 1.0 mm

451010

- ▶ Length: 130 mm
- ▶ Ø: 1.3 mm

451013



451010

URINARY CATHETER FOR MALE DOGS

Squeeze the water sachet before use to activate the hydrophilic coating.

- ▶ With Luer-Lock cone
- ▶ Extra long
- ▶ Soft and dimensionally stable
- ▶ Individually sterile packaged
- ▶ Pack of 5

- ▶ 6 Fr x 20"
- ▶ Length: 500 mm, Ø: 2 mm


452020

- ▶ 8 Fr x 20"
- ▶ Length: 500 mm, Ø: 2.6 mm

452026

- ▶ 10 Fr x 20"
- ▶ Length: 500 mm, Ø: 3.3 mm

452033



452020

BD PHASEAL™ ADAPTOR

ORDER
ONLINE

www.eickemeyer.com

- ▶ The BD PhaSeal™ adaptor is permanently connected to the drug vial. It creates and maintains a neutral pressure when air or liquid are injected into the patient.
- ▶ The only clinically proven closed-system for safe drug transfer
- ▶ Expansion chamber retains aerosols and vapours
- ▶ Double membrane provides dry and leakproof connections
- ▶ Needle Safe design allows complete drug removal from the bottle
- ▶ Includes a drug vial adaptor for closed reconstitution and pressure equalisation
- ▶ A component of the BD PhaSeal™ closed system transfer product range
- ▶ 1 Unit



409000

Adaptor P14 - Blue

- ▶ Neck size: 13 mm
- ▶ Equalizing capacity of 20 mL of air

409000

Adaptor P21 - Red

- ▶ Neck size: 20 mm
- ▶ Equalizing capacity of 60 mL of air

409001

Adaptor P28 - White

- ▶ Neck size: 28 mm
- ▶ Equalizing capacity of 60 mL of air

409004

Adaptor P50 - Green

- ▶ Neck size: 20 mm
- ▶ Equalizing capacity of 50 mL of air

409005

Adaptor P53 - Violet

- ▶ Neck size: 32 mm
- ▶ Equalizing capacity of 50 mL of air

409006

BD PHASEAL™ INFUSION

Part of the BD PhaSeal™ closed system transfer product range.

Infusion Adaptor C70

- ▶ IV connection with a built-in PhaSeal Connector for closed addition to IV bag
- ▶ Total length: 16 cm
- ▶ Package: 1 pc
- ▶ Carton: 50 pcs

409016



Infusion Adaptor C100

- ▶ Connection between IV bag and non vented infusion set
- ▶ Built-in PhaSeal Connector for closed addition to IV bag
- ▶ PVC free
- ▶ Package: 1 pc
- ▶ Carton: 200 pcs

409020



BD PhaSeal™ Infusion Set C50

- ▶ Standard IV infusion set with a built-in PhaSeal Connector for closed addition to IV bag
- ▶ Total length: 188 cm
- ▶ Priming volume: 17.5 ml
- ▶ Inline filter 15 µm. Drops/mL: 20
- ▶ Package: 1 pc
- ▶ Carton: 80 pcs

409009



BD PhaSeal™ Secondary Infusion Set C61

- ▶ Secondary set with a built-in PhaSeal Connector for closed addition to IV bag
- ▶ Total length: 49 cm
- ▶ Priming volume: 3.0 ml
- ▶ Package 1 pc
- ▶ Carton: 120 pcs

409011



BD PHASEAL™ CONNECTORS

Connector Luer Lock C35

- ▶ For I.V. Line access
- ▶ Female Luer-Lock connection

409003



Connector Luer Lock C45

Part of the BD PhaSeal™ closed system transfer product range.

- ▶ For I.V. Line access
- ▶ Female Luer-Lock connection

409010



M70 Connector Cap

Accessory for the BD PhaSeal™ closed drug transfer system. The connector cap (M70) was specifically designed for the infusion adapter (C70)

409022



Y-Site Connector C80

- ▶ Integrated connector with extension hose
- ▶ Length: 21 cm
- ▶ Volume: 1.6 ml

409017



M70 L Connector C90

Provides a dry connection for closed addition of hazardous drugs to semi-rigid plastic containers with rubber ports.

409007



BD PHASEAL™ INJECTORS

- ▶ Part of the BD PhaSeal™ closed system transfer product range
- ▶ Drug transfer device that attaches to a disposable syringe or IV tubing
- ▶ Drug transfer device that attaches to any disposable syringe or IV tubing with a Luer-Lock connection.
- ▶ Sold individually

BD PhaSeal™ Injector Luer Lock N35C

409002



BD PhaSeal™ Injector Luer Lock N30C

409008



BD PhaSeal™ Injector Luer Lock N35

409013



SURGICAL MASK

- ▶ Three ply
- ▶ Elastic band
- ▶ Pack of 50 pieces

050319



RESPIRATORY MASK FFP2

- ▶ Protecting both the wearer and the surrounding people
- ▶ Comfortable and easy to wear thanks to extra soft ear loops
- ▶ Adaptable (integrated) nose clip
- ▶ CE certified according to EN 149: 2001 + A1:2009
- ▶ Made in Germany
- ▶ Pack of 20 pieces

050324



ORDER
ONLINE

www.eickemeyer.com

FACE SHIELD

- ▶ Significantly decreases contamination risk
- ▶ Transparent polycarbonate with a wide angle of view
- ▶ Can be autoclaved at 121 °C for 20 minutes
- ▶ Forehead band consists of a highly-absorptive nonwoven fabric
- ▶ Pack of 10 pieces

940010



OVERSHOES

- ▶ Suitable for barrier nursing
- ▶ Disposable
- ▶ Pack of 100 pieces

050321



SEMPERMED CLASSIC SURGICAL GLOVE

- ▶ Sterile surgical glove with a high resistance to tearing and a reliable impermeability
- ▶ Natural rubber latex, hypo-allergenic
- ▶ Micro-rough surface ensures excellent sensitivity
- ▶ Sterile, in paper peel pouch
- ▶ Powdered
- ▶ 70 pairs per box

Size 6 Size 7 1/2
930760 930775

Size 6 1/2 Size 8
930765 930780

Size 7 Size 8 1/2
930770 930785



930770

STAINLESS STEEL GLOVE DISPENSER

With this wall-mountable glove dispenser, your examination gloves are always within easy reach. Even when the box is almost empty, the holder allows for easy and hygienic removal of the gloves.

- ▶ Fits 3 packs of nitrile examination gloves (item no. 050313 – 050317)
- ▶ Material: Stainless steel
- ▶ For wall mounting (incl. mounting hardware)
- ▶ Dimensions (in cm): W 21.2 x H 33.5 x D 6.8

050320



EXAMINATION GLOVE - LARGE ANIMAL

- ▶ Suitable for rectal examination in cows and horses
- ▶ High resistance to tearing
- ▶ Pack of 100 pieces

930191



NITRILE EXAMINATION GLOVES

- ▶ Latex free
- ▶ Powder free
- ▶ Ambidextrous
- ▶ Textured fingertips
- ▶ Non-sterile
- ▶ Colour: blue
- ▶ AQL: 1.5
- ▶ Thickness: palm 0.07 mm, finger 0.09 mm
- ▶ Total length: 240 mm
- ▶ Weight: 3.5 g
- ▶ Pack of 100 pieces

Size XS

050313

Size S

050314

Size M

050315

Size L

050316

Size XL

050317



BARRIER GOWN

ORDER
ONLINE

www.eickemeyer.com

- ▶ During procedures there is a high risk of being splashed by body fluids
- ▶ Our gowns are full length wrap-around with full length covering for arms but are soft and breathable
- ▶ They are fully disposable, which is ideal for preventing cross contamination between patients
- ▶ Excellent for protective barrier nursing procedures
- ▶ Disposable
- ▶ One size fits all

White

- ▶ Nonwoven
- ▶ Buttons on front

940752

Forest Green

- ▶ Nonwoven
- ▶ Arm cuffs and belt
- ▶ Pack of 10 pieces

940740

Green

- ▶ Plastic
- ▶ Adhesive closure

940275



940740

EICKWRAP COHESIVE BANDAGES

EickWrap bandage's self-adhesive surface means that they adhere to each other well but do not fix to skin, hair or other textiles. EickWrap does not constrict blood circulation or joint movements. The bandages are best suited for horse legs and hooves or dog and cat paws.

- ▶ Highest quality cohesive bandage
- ▶ Conforms itself to the shape of the area where it is applied – helping to keep dressing in place
- ▶ Self-adhesive but, it will not stick to the skin
- ▶ Contains latex
- ▶ Allows air to circulate over the skin and around the dressing
- ▶ Individually packaged – assorted colours
- ▶ Length: 4.5 m/roll

- ▶ Width: 50 mm
- ▶ Pack of 36 rolls

050750

- ▶ Width: 75 mm
- ▶ Pack of 24 rolls

050751

- ▶ Width: 100 mm
- ▶ Pack of 18 rolls

050752



050750

050751

EICKWRAP FUN COHESIVE BANDAGES

The tried and tested EickWrap bandages are now also available with fun, colourful motifs. The bandages have a self-adhesive surface structure so that they stick to one another without any problems.

- ▶ 6 Colourful Motifs: Choose from a variety of eye-catching designs
- ▶ Highest Quality Cohesive Bandage: Reliable and durable
- ▶ Conforms to Shape: Adapts to the area where it is applied
- ▶ Breathable: Allows air circulation over the skin and around the dressing
- ▶ Individually Packaged: Convenient and hygienic
- ▶ Contains latex
- ▶ Length: 4.5m/roll
- ▶ Width: 5 cm
- ▶ Box of 36 pieces

050758



EICKPAD UNDERCAST PADDING

- ▶ Breathable and temperature regulating
- ▶ Areas of use: protection of exposed bone or nerve points underneath bandages and padding of splints and casts, in combination with wound dressing as breathable padding for wound management
- ▶ Material: 100 % Polyester

- ▶ 6 cm x 2.7 m
- ▶ Pack of 12 pieces

050755

- ▶ 10 cm x 2.7 m
- ▶ Pack of 12 pieces

050756

- ▶ 15 cm x 2.7 m
- ▶ Pack of 6 pieces

050757



050756

GUTTER SPLINT

- ▶ For temporary immobilisation of fractures
- ▶ Strong and lightweight plastic, green
- ▶ Comfortable, non-irritating foam padding
- ▶ Anatomically shaped, fully radiolucent
- ▶ Scoop end to be placed at toes

Dimensions (in cm):
L 20 x W 3

951720

Dimensions (in cm):
L 40 x W 5

951740

Dimensions (in cm):
L 30 x W 4

951730



951740

951730

951720

EICKSPLINT

- ▶ Malleable aluminium splint that can be formed into any desired shape and can be easily cut to the desired size
- ▶ Sturdy aluminium-rail core with foam padding
- ▶ Dimensions (in cm): L 100 x W 11

951811



ABDOMINAL BANDAGE

ORDER
ONLINE

www.eickemeyer.com

- ▶ Soft and non abrasive on the skin
- ▶ 100 % cotton
- ▶ Washable at 95 °C
- ▶ Adjustable Velcro strap to ensure a perfect fit
- ▶ Tailored fit with 7 different sizes
- ▶ Original only from EICKEMEYER®

Size 1

- ▶ Suitable for cats
- ▶ Length: 32 cm

950001

Size 2

- ▶ Suitable for Poodles
- ▶ Length: 38 cm

950002

Size 3

- ▶ Suitable for Setters
- ▶ Length: 47 cm

950003

Size 4

- ▶ Suitable for Labradors
- ▶ Length: 51 cm

950004

Size 5

- ▶ Suitable for German Shepherds
- ▶ Length: 60 cm

950005

Size 6

- ▶ Suitable for Great Danes
- ▶ Length: 67 cm
- ▶ Width: 25 cm

950006

Size 7

- ▶ Suitable for St. Bernards
- ▶ Length: 67 cm
- ▶ Width: 25 cm

950007



application

GAUZE SWABS

- ▶ 100 % cotton
- ▶ Non-sterile
- ▶ 8 ply
- ▶ 10 packs of 100 each (packed in a box)

Dimensions (in cm):

L 5 x W 5

050080

Dimensions (in cm):

L 10 x W 10

050082

Dimensions (in cm):

L 7.5 x W 7.5

050081



050080

EICKEMEYER® PAPER TOWELS

- ▶ Cleaning tissue made from unbleached pulp
- ▶ Dimensions (in cm): W 40 x L 60
- ▶ 15 kg pack

953815



VETMEDCARE DOG BODY BASIC

The high-quality animal bodies from VetMedCare® are now available at EICKEMEYER®. They are used after operations and protect the wound from dirt, licking and scratching.

- ▶ Skin-friendly and easy-care material (90 % polyamide, 10 % elastane)
- ▶ Protects the skin against allergies and irritations
- ▶ Comfortable and hygienic as washable up to 60 °C
- ▶ Easy urination (extra opening for males) and defecation
- ▶ Developed in collaboration with veterinary surgeons
- ▶ Special flat and stretchy seams prevent pressure points
- ▶ Completely smooth abdominal area
- ▶ Breathable and extremely quick drying
- ▶ Patient comfort due to the enormous elasticity and ergonomic adjustment




VMC10018

VetMedCare Dog Body		
Item No.	Description	Size
VMC10008	Dog Body	XXXS unisex
VMC10009	Dog Body	XXS unisex
VMC10011	Dog Body	XS bitch
VMC10021	Dog Body	XS male
VMC10017	Dog Body	XS plus bitch
VMC10027	Dog Body	XS plus male
VMC10012	Dog Body	S bitch
VMC10022	Dog Body	S male
VMC10018	Dog Body	S plus bitch
VMC10028	Dog Body	S plus male
VMC10013	Dog Body	M bitch
VMC10023	Dog Body	M male
VMC10014	Dog Body	L bitch
VMC10024	Dog Body	L male
VMC10015	Dog Body	XL bitch
VMC10025	Dog Body	XL male
VMC10016	Dog Body	XXL bitch
VMC10026	Dog Body	XXL male

VETMEDCARE SIZE CHART BODY

- ▶ The body should be worn very tightly, this is particularly important for bitches so that the body does not get dirty when urinating
- ▶ Breed information is to be considered as a guide; If you are hesitating between two sizes due to the length, it is advisable to choose the larger size for corpulent animals, and the smaller size can be used for slim animals



Size Chart Body		
Size	Breed (example)	Length (in cm)
XXXS	Baby Kitten, Guinea Pig, Rat	24 – 26
XXS	Small Cat, Shih Tzu	27 – 31
XS	Chihuahua, Lhasa Apso, Small Yorkshire Terrier	32 – 38
XS plus	Yorkshire Terrier, Lhasa Apso, Papillon, Longhair Chihuahua	39 – 44
S	Maltese, Cavalier King Charles, Miniature Dachshund	45 – 50
S plus	French Bulldog, Pug, Jack Russel, Dachshund, Japanese Spitz	51 – 54
M	Schnauzer, Shorthair Border Collie, Cocker Spaniel	55 – 62
L	American Bulldog, Golden Retriever, Labrador, Doberman	63 – 74
XL	German Shepherd, Rottweiler, Bernese Mountain Dog	75 – 79
XXL	Saint Bernard, Mastiff, Newfoundland, Great Dane	from 80

VETMEDCARE DOG BODY WITH LEGS

The dog body with legs has 4 extra long legs already sewn on, which can be individually adjusted in length due to the run-proof material.

- ▶ Optimal wound protection for leg injuries, e.g. after cruciate ligament surgery
- ▶ Also ideal after castration, as the legs provide a better closure
- ▶ Can also be used for paw injuries, the overly long legs can be turned inwards
- ▶ Tick protection without chemicals
- ▶ Comfortable wound protection (leak protection) without collar
- ▶ Comfortable and hygienic, as washable up to 60 °C



VMC12018



VetMedCare Dog Body with legs		
Item No.	Description	Size
VMC12008	Dog Body with legs	XXXS unisex
VMC12009	Dog Body with legs	XXS unisex
VMC12011	Dog Body with legs	XS bitch
VMC12021	Dog Body with legs	XS male
VMC12017	Dog Body with legs	XS plus bitch
VMC12027	Dog Body with legs	XS plus male
VMC12012	Dog Body with legs	S bitch
VMC12022	Dog Body with legs	S male
VMC12018	Dog Body with legs	S plus bitch
VMC12028	Dog Body with legs	S plus male
VMC12013	Dog Body with legs	M bitch
VMC12023	Dog Body with legs	M male
VMC12014	Dog Body with legs	L bitch
VMC12024	Dog Body with legs	L male
VMC12015	Dog Body with legs	XL bitch
VMC12025	Dog Body with legs	XL male
VMC12016	Dog Body with legs	XXL bitch
VMC12026	Dog Body with legs	XXL male

VETMEDCARE DOG BODY WITH LEGS AND ZIP – UNISEX

With zip for quick and easy dressing and undressing. Thanks to the run-proof material, the additional opening for the male dog can be easily cut with scissors and the leg length can be individually adjusted.

- ▶ With loop to fix the collar



VMC12053

VetMedCare Dog Body with legs and zip – unisex		
Item No.	Description	Size
VMC12048	Dog Body with legs and zip	XXXS
VMC12049	Dog Body with legs and zip	XXS
VMC12051	Dog Body with legs and zip	XS
VMC12057	Dog Body with legs and zip	XS plus
VMC12052	Dog Body with legs and zip	S
VMC12058	Dog Body with legs and zip	S plus
VMC12053	Dog Body with legs and zip	M
VMC12054	Dog Body with legs and zip	L
VMC12055	Dog Body with legs and zip	XL
VMC12056	Dog Body with legs and zip	XXL

VETMEDCARE CAT BODY – UNISEX

The ideal solution for all cats.

- ▶ Skin-friendly and easy-care material (90 % polyamide, 10 % elastane)
- ▶ Protects the skin against allergies and irritations
- ▶ Comfortable and hygienic as washable up to 60 °C
- ▶ Highly elastic material, is worn close-fitting and is therefore like a second skin
- ▶ Developed in collaboration with Veterinary surgeons
- ▶ Special flat and stretchy seams prevent unnecessary pressure points
- ▶ Completely smooth abdominal area
- ▶ Breathable and extremely quick drying
- ▶ Patient comfort due to the enormous elasticity and ergonomic adjustment
- ▶ With extra long neck
- ▶ With 4 short legs
- ▶ Dual label loops for collar attachment
- ▶ Comfortable fit, easy to put on and take off

VetMedCare Cat Body – unisex		
Item No.	Description	Size
VMC14011	Cat Body	XS
VMC14012	Cat Body	S
VMC14017	Cat Body	XS plus



VMC14012

VETMEDCARE SAFETY TUBE – UNISEX

The 1 m long, individually cut tubular fabric (run-proof material) is the ideal wound protection for small animals.

Simply cut the tube to the desired length and cut out the required openings – the body for your ferret, guinea pig, rabbit, hamster, koi, etc. is complete.

Attached to the dog or cat body, the safety tube can also be used as rod protection for use over the tail.

- ▶ Skin-friendly and easy-care material (90 % polyamide, 10 % elastane)
- ▶ Comfortable and hygienic, as washable up to 60° C

VetMedCare Safety Tubes – unisex		
Item No.	Description	Size
VMC21011	Safety Tube	S (3,5 cm)
VMC21012	Safety Tube	M (5 cm)
VMC21013	Safety Tube	L (8 cm)
VMC21014	Safety Tube	Conical (top 13.5 cm, bottom 4 cm)
VMC21010	Safety Tube	Set of 4, S + M + L + Conical



VMC21013

VETMEDCARE SAFETY PANTS – UNISEX

**ORDER
ONLINE**

www.eickemeyer.com

Suitable protection for a wide variety of ailments, as a stand alone product or in conjunction with the Cat and Dog Body range.

- ▶ Ideal for castration
- ▶ In cases of heat or incontinence, sanitary/incontinence pads can be inserted under the body
- ▶ With cord to attach to chest harness
- ▶ Skin-friendly and easy-care material (90 % polyamide, 10 % elastane)
- ▶ Comfortable and hygienic, washable up to 60 °C



VetMedCare Safety Pants – unisex		
Item No.	Description	Size
VMC13011	Safety Pants	XS / S
VMC13012	Safety Pants	S plus
VMC13013	Safety Pants	M / L
VMC13015	Safety Pants	XL / XXL

Size chart Safety Pants – unisex	
Size	Breed (example)
XS / S	Chihuahua, Papillon, Small Yorkshire Terrier, Lhasa Apso, Cavalier King Charles, Miniature Dachshund
S plus	French Bulldog, Pug, Maltese, Jack Russel, Japanese Spitz
M / L	Schnauzer, Border Collie, Cocker Spaniel, Labrador, Golden Retriever, American Bulldog
XL / XXL	Rottweiler, Bernese Mountain Dog, St. Bernard, Mastiff, Newfoundland



VMC13013

VETMEDCARE PRESS STUD

Press studs for attaching a safety tube (leg) to the dog or cat body.

- ▶ Colour: black

VetMedCare Press Studs		
Item No.	Description	Size
VMC41095	Press Studs, 5 pcs	21 mm
VMC41096	Press Studs, 6 pcs	11 mm



VMC41095

PATIENT SUPPORT FOR ANIMALS

- ▶ Soft top made of non-woven
- ▶ Material Base: Polyethylene
- ▶ Very absorbent (suction capacity 2,000 ml)
- ▶ Pack of 30 pieces
- ▶ Dimensions (in cm): L 90 x W 60

050070



SCAVENGING SYSTEM EVACUATION TUBE

- ▶ Ø 22 mm
- ▶ Supplied per metre
- ▶ Transparent

212729



EICKSORBER ANAESTHETIC GAS ABSORBER

- ▶ Activated carbon filter to absorb organic anaesthetic gases.
- ▶ Connects to pressure relief valve. Thus preventing unwanted gases from entering room atmosphere
- ▶ Replace filters after exceeding 1,000 g
- ▶ Includes 2 m exhaust hose (Ø 22 mm)

213488

Single Unit

- ▶ Without exhaust hose

21348801

Pack of 6 pieces

- ▶ Includes 2 m exhaust hose (Ø 22 mm)

21348806



SODA LIME ANAESTHESIA

- ▶ Medical grade Soda Lime is used to remove carbon dioxide in closed anaesthesia circuits
- ▶ Use with rebreathing circuits to prevent rebreathing of CO₂
- ▶ Low dust particles
- ▶ Colour indicator: white to violet
- ▶ Content: 5 l

215206



ANAESTHETIC CIRCUIT TUBING

- ▶ Y-Connector
- ▶ Transparent
- ▶ Smoothbore
- ▶ Ø 22 mm
- ▶ Length: 1.2 m


212718



AYRE'S T-PIECE BREATHING CIRCUIT

- ▶ For patients < 7 kg
- ▶ Pressure relief valve
- ▶ Length: 1.8 m
- ▶ Includes: 0.5 l rebreathing bag

21309514





EICKTUBE PVC

**ORDER
ONLINE**

www.eickemeyer.com

- ▶ Valve for Luer and Luer lock fittings
- ▶ Flattened, atraumatic tip
- ▶ Includes connector

EickTube Starter Set				
Item No.	Description			
215800	Consisting of 23 tracheal tubes up to inner diameter of 10.0			
Single				
Item No.	Quantity	Description	Ø inner mm	Ø outer mm
215920	25 pcs	Cuffless	2.0	2.9
215820	25 pcs	Cuff (Ø 9 mm)	2.0	2.9
215925	25 pcs	Cuffless	2.5	3.5
215825	25 pcs	Cuff (Ø 9 mm)	2.5	3.5
215930	25 pcs	Cuffless	3.0	4.2
215830	25 pcs	Cuff (Ø 9 mm)	3.0	4.2
215935	25 pcs	Cuffless	3.5	4.9
215835	25 pcs	Cuff (Ø 11 mm)	3.5	4.9
215940	25 pcs	Cuffless	4.0	5.5
215840	25 pcs	Cuff (Ø 13 mm)	4.0	5.5
215945	25 pcs	Cuffless	4.5	6.2
215845	25 pcs	Cuff (Ø 13 mm)	4.5	6.2
215850	25 pcs	Cuff (Ø 13 mm)	5.0	6.9
215855	25 pcs	Cuff (Ø 18 mm)	5.5	7.5
215860	25 pcs	Cuff (Ø 22 mm)	6.0	8.2
215865	25 pcs	Cuff (Ø 22 mm)	6.5	8.7
215870	25 pcs	Cuff (Ø 25 mm)	7.0	9.3
215875	25 pcs	Cuff (Ø 25 mm)	7.5	10.0
215880	25 pcs	Cuff (Ø 27 mm)	8.0	10.7
215885	25 pcs	Cuff (Ø 27 mm)	8.5	11.3
215890	25 pcs	Cuff (Ø 29 mm)	9.0	12.0
215895	25 pcs	Cuff (Ø 29 mm)	9.5	12.7
215910	10 pcs	Cuff (Ø 29 mm)	10.0	13.3
215911	10 pcs	Cuff (Ø 35 mm)	11.0	14.7
215912	10 pcs	Cuff (Ø 52 mm)	12.0	16.0
215913	10 pcs	Cuff (Ø 52 mm)	13.0	17.3



215845

ENDOTRACHEAL TUBE LUBRICANT SPRAY

Silicone based lubricant that protects tubes.

- ▶ Used for maintenance of tracheal tubes, probes and catheters
- ▶ Used when storing tracheal tubes and to prolong storage lifespan
- ▶ WARNING: Not to be used as a lubricant or to ease insertion of tubes or catheters
- ▶ It is CFC-free and only suitable for lubricating moving parts and preventing rusting of metal surgical instruments

214600



MADGIC™ DISPOSABLE DRUG DISPENSER

- ▶ MADGIC™ is a form-stable tool for the application of local anaesthetics
- ▶ Targeted application to the nasal / oropharynx and hypopharynx mucosa
- ▶ The Disposable Drug Dispenser consists of a flexible and at the same time stable cannula, which atomizes drugs into a fine mist of 30–100 micron particles

Single Unit

05040001



EICKFIL NEEDLE-SUTURE COMBINATIONS AND SPECIAL SUTURES WITHOUT NEEDLE

Round bodied

- ▶ Braided and Coated Absorbable Suture Material
- ▶ Coating: Calcium Stearate & Polycaprolactone
- ▶ Material: Coated Polyglycolic Acid
- ▶ Absorption by hydrolysis
- ▶ Tensile strength: 84 % at 14 days, 23 % at 28 days
- ▶ Total absorption: 60 to 90 days
- ▶ Colour: Violet
- ▶ Taper Point: 1/2 Circle
- ▶ Pack of 12 pieces
- ▶ Thread length: 70 cm

Item No.	mm	USP	Metric
050607	26	3/0	2
050615	30	3/0	2
050608	26	2/0	3
050616	30	2/0	3
050617	30	0	3.5
050618	30	1	4

Reverse cutting

- ▶ Braided and Coated Absorbable Suture Material
- ▶ Material: Coated Polyglycolic Acid
- ▶ Coating: Calcium Stearate & Polycaprolactone
- ▶ Absorption by hydrolysis
- ▶ Tensile strength: 84 % at 14 days, 23 % at 28 days
- ▶ Total absorption: 60 to 90 days
- ▶ Reverse Cutting: 3/8 Circle
- ▶ Colour: Violet
- ▶ Pack of 12 pieces
- ▶ Thread length: 70 cm

Item No.	mm	USP	Metric
050650	19	5/0	1
050651	19	4/0	1.5
050658	24	4/0	1.5
050652	19	3/0	2
050659	24	3/0	2
050665	30	3/0	2
050653	19	2/0	3
050660	24	2/0	3
050666	30	2/0	3
050661	24	0	3.5
050667	30	0	3.5
050668	30	1	4

Special Sutures without Needle

- ▶ Material: 100 % polyglycolic acid – Braided & Coated
- ▶ Colour: Violet
- ▶ Coating: Calcium Stearate & Polycaprolactone
- ▶ Sterilization: Ethylene Oxide
- ▶ Shelf life: 5 years
- ▶ Thread length: 2.5 m

Item No.	Metric	USP
051072	2	3/0
051073	3	2/0
051074	3/4	0
051075	4	1
051076	5	2

Features:

- ▶ Absorption by hydrolysis, within 60–90 days
- ▶ Minimal tissue reaction
- ▶ Coating ensures smooth tissue passage
- ▶ Outstanding pliability
- ▶ Individually packaged – sterile
- ▶ Pack of 12 pieces



051075

EICKSKIN SKIN STAPLER

**ORDER
ONLINE**

www.eickemeyer.com

EickSkin Skin Stapler

The EickSkin Disposable Skin Stapler is a simple and affordable stapler that provides consistent and reliable performance. Designed to fit comfortably in your hand, it enhances control of skin closure reducing surgery time as well as minimising infections and scarring. In addition the EickSkin Staples can be easily removed causing no pain to the patient.

- ▶ Complete with 35 staples (7 mm)
- ▶ Individually packaged – sterile
- ▶ Pack of 10 pieces

157312



157312

Skin Staple Remover

- ▶ Dedicated surgical staple remover

157309

DISPOSABLE SKIN PUNCHES

- ▶ Suitable for skin and tissue biopsy
- ▶ Disposable
- ▶ Individually packaged – sterile
- ▶ Pack of 10 pieces

Ø 2 mm

179320

Ø 6 mm

179360

Ø 4 mm

179340

Ø 8 mm

179380



179340

EICKEMEYER® SCALPEL BLADES

- ▶ Fits STANDARD scalpel handles No. 3 (Figs. 10–15) and No. 4 (Figs. 18–23)
- ▶ Superior cutting quality thanks to precision grinding
- ▶ Increased safety through precise cutting
- ▶ Wide selection of blades
- ▶ Individually sterile-packed
- ▶ Bulk pack of 100

Fig. 10

100810

Fig. 11

100811

Fig. 15

100815

Fig. 20

100820

Fig. 21

100821

Fig. 22

100822

Fig. 23

100823



100823

MICRO SCALPEL BLADES

- ▶ For ophthalmic, neurosurgical, and microsurgical use
- ▶ Extra-strong blades for precise cutting
- ▶ Extremely sharp
- ▶ Pack of 25

Fig. 63

101863

Fig. 64

101864

Fig. 65

101865

Fig. 69

101869



100823

EICKFIL OPHTHALMIC SUTURE MATERIAL

PGA

- ▶ Type: Braided and coated multifilament
- ▶ Coating: Polycaprolactone and calcium stearate (1 %)
- ▶ Colour: Violet / Beige – undyed
- ▶ Tissue reaction: Minimal
- ▶ Absorption: The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days. Approximately 50 % of tensile strength remain after 21 days.
- ▶ Presentation: Box with 12 sutures

Characteristics

- ▶ High tensile strength
- ▶ Synthetic
- ▶ Breaks down by hydrolysis
- ▶ Braided structure
- ▶ Coated multifilament
- ▶ Hermetically sealed packing

- ▶ USP 9/0, EP 0.3
- ▶ 3/8 circle, DSP 6.4 mm, 200 µm
- ▶ Side cutting spatula
- ▶ Colour: violet
- ▶ Length: 30 cm

050889

- ▶ USP 3/0, EP 2
- ▶ 3/8 circle, DS 24 mm
- ▶ Reverse cutting
- ▶ Colour: undyed
- ▶ Length: 75 cm

050892

- ▶ USP 9/0, EP 0.3
- ▶ 3/8 circle, 2 x DSP 6.4 mm, 200 µm
- ▶ Side cutting spatula double armed
- ▶ Colour: violet
- ▶ Length: 30 cm

050895

- ▶ USP 6/0, EP 0.7
- ▶ 1/4 circle, 2 x VSP 8 mm, 300 µm
- ▶ Side cutting spatula double armed
- ▶ Colour: violet
- ▶ Length: 45 cm

050898

- ▶ USP 6/0, EP 0.7
- ▶ 3/8 circle, DS 12 mm
- ▶ Reverse cutting
- ▶ Colour: violet
- ▶ Length: 45 cm

050890

- ▶ USP 6/0, EP 0.7
- ▶ 1/2 circle, HR 17 mm
- ▶ Round bodied taper point
- ▶ Colour: violet
- ▶ Length: 75 cm

050893

- ▶ USP 8/0, EP 0.4
- ▶ 3/8 circle, 2 x DSP 6 mm, 200 µm
- ▶ Side cutting spatula double armed
- ▶ Colour: violet
- ▶ Length: 45 cm

050896

- ▶ USP 3/0, EP 2
- ▶ 3/8 circle, DS 19 mm easy pass
- ▶ Reverse cutting
- ▶ Colour: undyed
- ▶ Length: 50 cm

050891

- ▶ USP 10/0, EP 0.2
- ▶ 3/8 circle, 2 x DSP 6.2 mm, 150 µm
- ▶ Side cutting spatula double armed
- ▶ Colour: violet
- ▶ Length: 30 cm

050894

- ▶ USP 7/0, EP 0.5
- ▶ 3/8 circle, 2 x DSP 8 mm, 200 µm
- ▶ Side cutting spatula double armed
- ▶ Colour: violet
- ▶ Length: 45 cm

050897



Side cutting spatula
3/8 circle



Reverse cutting
3/8 circle



Round bodied taper point
1/2 circle



Side cutting spatula double armed
3/8 circle



Side cutting spatula double armed
1/4 circle

Nylon

- ▶ Type: Monofilament
- ▶ Composition: Extrusion of polyamide 6.0 or 6.6
- ▶ Coating: None
- ▶ Colour: Black
- ▶ Tissue reaction: Minimal
- ▶ Absorption: Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10 % a year by rupture of chemical links (hydrolytic action).
- ▶ Presentation: Box with 12 sutures

Characteristics

- ▶ High tensile strength, non absorbable
- ▶ Synthetic
- ▶ Monofilament with microscopic uniform diameter
- ▶ Pliability and softness
- ▶ Hermetically sealed packing

- | | |
|--------------------------------------|--------------------------------------|
| ▶ USP 10/0, EP 0.2 | ▶ USP 9/0, EP 0.3 |
| ▶ 3/8 circle, 2 x DSP 6.4 mm, 220 µm | ▶ 3/8 circle, 2 x DSP 6.2 mm, 150 µm |
| ▶ Side cutting spatula double armed | ▶ Side cutting spatula double armed |
| ▶ Monofilament | ▶ Monofilament |
| ▶ Colour: black | ▶ Colour: black |
| ▶ Length: 45 cm | ▶ Length: 45 cm |
| ▶ Quantity: 12 sutures | ▶ Quantity: 12 sutures |

050899

050900

Silk

- ▶ Type: Braided multifilament
- ▶ Composition: Braided fibres from the cocoon of the silkworm
- ▶ Coating: Waxes silicone
- ▶ Colour: Black
- ▶ Tissue reaction: Moderate
- ▶ Absorption: Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.
- ▶ Presentation: Box with 12 sutures

Characteristics

- ▶ High tensile strength, non absorbable
- ▶ Braided or twisted structure
- ▶ Coated multifilament
- ▶ Hermetically sealed packing

- ▶ USP 4/0, EP 1.5
- ▶ 3/8 circle, DS 19 mm
- ▶ Reverse cutting
- ▶ Braided multifilament
- ▶ Colour: black
- ▶ Length: 45 cm
- ▶ Quantity: 12 sutures

050901



050889

BAUSCH & LOMB PUREVISION® 2 BANDAGE CONTACT LENSES

Contact bandage lenses act as a protective shield on the cornea. They permit corneal repair and improve patient comfort by preventing contact between damaged cornea and eyelids.

- ▶ High definition optics offers a clear and sharp view
Excellent stabilisation
- ▶ Comfortable contact lens wear for up to 30 days and nights
- ▶ Ensures an optimal supply of oxygen to the cornea and makes eyes look clear
- ▶ Packaging: 6 lenses per box

173540




SCHIRMER TEAR TEST STRIPS

SCHIRMER Tear test can be used to measure basal and reflexive tear production in non-anaesthetised eye.

- ▶ Blue colour mark tear front indicator
- ▶ Physiological range in dogs: 15 – 25 mm
- ▶ Individually packed – sterile
- ▶ Pack of 100 pieces

173592



CAUTERY PEN

- ▶ Precise cauterisation thanks to the very fine tip (28 mm)
- ▶ Easy to use
- ▶ Operating temperature: 600 °C
- ▶ Individually packaged – sterile
- ▶ For single use only

173430



SUGI SWABS

- ▶ Sugi absorbent materials are especially suitable for dabbing those areas where special care must be taken (e. g. neural structures and ophthalmic tissue)
- ▶ Sugi absorbent materials are manufactured from cellulose material
- ▶ The material undergoes soft / elastic expansion immediately upon contacting moisture
- ▶ The absorptive capacity is approximately 17 times its own weight
- ▶ Length: 70 mm

Wedge Tip

- ▶ 7 mm
- ▶ Pack of 200 pieces
(each 10 pieces sterile packed)

166520

Cylindrical Shape

- ▶ Pack of 500 pieces

166530



OCU-FILM® TONO-PEN TIP COVERS

- ▶ Minimises the risk of cross-contamination
- ▶ Pack of 150 pieces

173565



TONOVET PROBES

- ▶ For TONOVET (1736005) and TONOVET Plus (1736007)
- ▶ Pack of 100 pieces

1736006



**ORDER
ONLINE**
www.eickemeyer.com

AUTOCLAVE STERILISATION REEL

- ▶ Tear resistant
- ▶ Bacterial resistant
- ▶ Transparent on one side
- ▶ Heat sealable with indicator
- ▶ Reel length: 200 metres

Width: 50 mm

502005

Width: 150 mm

502015

Width: 70 mm

502007

Width: 200 mm

502020

Width: 100 mm

502010

Width: 250 mm

502025



502015

AUTOCLAVE STERILISATION POUCH

- ▶ Self-sealing with indicator
- ▶ Pack of 200 pieces

Dimensions (in mm):
L 250 x W 90

502209

Dimensions (in mm):
L 260 x W 140

502213



502209

INSTRUMENT LUBRICANT

- ▶ For maintenance of instruments prior to sterilisation
- ▶ Content: 500 ml

563707



DISPOSABLE POLISHING HEADS (HP)

Optimal polishing thanks to an oscillating movement. This also prevents excessive heat from developing on the tooth and the deformation of whisker hair.

- ▶ For straight handpiece
- ▶ Max. speed 3,000 rpm
- ▶ For single use
- ▶ Pack of 100 pieces

174066



POLISHING CUPS (RA)

Perfect shiny teeth with a highly polished enamel surface.

- ▶ Compatible with Contra Angle Handpiece
- ▶ Pack of 5 pieces

176787



EICKEMEYER® PROPHY PASTE

- ▶ For polishing teeth after removing tartar
- ▶ Medium grain
- ▶ Peppermint flavour
- ▶ Content: 200 g

176782



DENTAL HANDPIECE ANTI DUST SPRAY

- ▶ Suitable to remove dust from handpieces
- ▶ CFC free
- ▶ Includes spray head
- ▶ Content: 400 ml

174880



DENTAL HANDPIECE LUBRICANT OIL

- ▶ A specially formulated, highly effective general handpiece and air motor lubricant
- ▶ Frequent application will assure hand piece and air motor performance and durability
- ▶ Suitable for high & low speed handpieces
- ▶ Content: 500 ml

174010



AIR MOTOR / TURBINE ADAPTER AND LOW SPEED E-TYPE ADAPTER

Low Speed E-Type Adapter

- ▶ Suitable for Dental Handpiece Lubricant Oil (174010)

17401155

Air Motor / Turbine Adapter

- ▶ For high speed handpieces

17401150



17401150



17401155

EICKCOLL DOG COLLARS

**ORDER
ONLINE**

www.eickemeyer.com

- ▶ Prevent animals from reaching wounds
- ▶ Easy to fit on patient
- ▶ Made of flexible polyethylene
- ▶ Supple and easy to handle

- ▶ Ø neck: 7 to 8 cm
- ▶ Length: 7.5 cm

220020

- ▶ Ø neck: 10 to 12 cm
- ▶ Length: 15 cm

220023

- ▶ Ø neck: 15 to 18 cm
- ▶ Length: 30 cm

220026

- ▶ Ø neck: 7 to 8 cm
- ▶ Length: 10 cm

220021

- ▶ Ø neck: 12 to 14 cm
- ▶ Length: 20 cm

220024

- ▶ Ø neck: 9 to 10.5 cm
- ▶ Length: 12 cm

220022

- ▶ Ø neck: 14 to 16 cm
- ▶ Length: 25 cm

220025



220025

EICKEMEYER® TICK REMOVER

- ▶ Suitable for pet owners
- ▶ Fast and safe removal of ticks
- ▶ Length: 70 mm

907321



SMALL ANIMAL PILL GIVER

- ▶ A convenient way to administer pills to small animals
- ▶ Pack of 12 pieces

951802



EICKEMEYER® ULTRASOUND GEL

The EICKEMEYER® Ultrasound Gel is a proven contact agent for ultrasound diagnostics and therapy.

- ▶ Free of bubbles and acoustically precise scans
- ▶ Free of Formaldehyde, perfume additives, dyes and UV stabilizers therefore minimal allergenic potency
- ▶ Not suitable for ophthalmological examinations

303936



EICKEMEYER® CLIPPER SPRAY

- ▶ Clipper blades require cleaning in-between each patient to provide the longest possible working life
- ▶ Dual action of disinfection and lubrication
- ▶ Content: 400 ml

902513



ENDOTRACHEAL TUBE LUBRICANT SPRAY

Silicone based lubricant that protects tubes.

- ▶ Used for maintenance of tracheal tubes, probes and catheters
- ▶ Used when storing tracheal tubes and to prolong storage lifespan
- ▶ **WARNING:** Not to be used as a lubricant or to ease insertion of tubes or catheters
- ▶ It is CFC-free and only suitable for lubricating moving parts and preventing rusting of metal surgical instruments

214600



SCHÜLKE ANTIFECT®

Schülke antifect® N liquid surface disinfectant is a rapid alcohol-based disinfection of non-invasive medical devices and surfaces.

- ▶ Very broad-spectrum effectiveness within a very short time
- ▶ Bactericidal – 1 min.
- ▶ Yeasticidal – 1 min.
- ▶ Limited virucidal – 30 sec.
- ▶ Leaves no streaks
- ▶ Ready to use
- ▶ 1 l bottle

307200



SCHÜLKE SPRAY PUMP

- ▶ Spray pump for liquids
- ▶ Fits 500 ml and 1 l bottles

307201



SCHÜLKE ANTIFECT®

Schülke antifect® N liquid surface disinfectant is a rapid alcohol-based disinfection of non-invasive medical devices and surfaces.

- ▶ Very broad-spectrum effectiveness within a very short time
- ▶ Bactericidal – 1 min.
- ▶ Yeasticidal – 1 min.
- ▶ Limited virucidal – 30 sec.
- ▶ Leaves no streaks
- ▶ Ready to use
- ▶ 5 l container

307205



Concentrate for basic disinfection and cleaning of medical devices and surfaces.

- ▶ Excellent cleaning performance
- ▶ Free from dyes and fragrances
- ▶ Complies with current EN standards
- ▶ Suitable for basic disinfection of surfaces (PT2) and medical devices
- ▶ Bactericidal: 0.5% (5 ml/l) – 30 min., 1% (10 ml/l) – 5 min.
- ▶ Yeastocidal: 0.5% (5 ml/l) – 30 min., 1% (10 ml/l) – 5 min.
- ▶ Limited virucidal activity: 0.5% (5 ml/l) – 30 min., 1% (10 ml/l) – 5 min.

Schülke Puresept® Surface Disinfectant

- ▶ 2 l bottle

307220



Schülke Puresept® Surface Disinfectant

- ▶ 5 l canister

307225



SCHÜLKE DESDERMAN®

Dye- and fragrance-free alcoholic hand rub for hygienic and surgical hand disinfection with nourishing panthenol and vitamin E. The integrated dosing pump allows for hygienic, precise, and drip-free dispensing of the product, ideal for use in high-traffic areas.

- ▶ Year-round product: reliably protects against infections all year round
- ▶ Broad spectrum of activity, including virucidal activity according to EN 14476
- ▶ Skin care formula with ProPanthenol Complex + Vitamin E
- ▶ Free from residual active ingredients
- ▶ Biocidal product
- ▶ Bactericidal – 15 sec.
- ▶ Mycobactericidal – 15 sec.
- ▶ Yeastocidal – 15 sec.
- ▶ Fungicidal – 60 sec.
- ▶ Limited virucidal – 15 sec.
- ▶ Virucidal – 30 sec.

Schülke desderman® Care Hand Disinfectant

- ▶ 1 l bottle

307230



Schülke Desderman® Care Hand Sanitizer bottle with dosing pump

- 1 l bottle

307231



Schülke Desderman® care Hand disinfectant

- 5 l canister

307235



SCHÜLKE GIGAZYME®

Schülke gigazyme® Instrument cleaner is an enzymatic cleaner for manual reprocessing of medical devices, such as flexible endoscopes and surgical instruments.

- ▶ First-class cleaning power thanks to three enzymes
- ▶ Excellent material compatibility thanks to a neutral pH value
- ▶ Low foaming thanks to non-ionic surfactants
- ▶ Optimal combination of high-quality enzymes (protease = breaks down proteins, lipase = dissolves fatty soils, amylase = removes polysaccharides)
- ▶ Excellent cleaning of all types of thermostable and thermolabile medical devices with excellent material compatibility
- ▶ 2 l bottle

307300



Schülke gigazyme® Instrument cleaner is an enzymatic cleaner for manual reprocessing of medical devices, such as flexible endoscopes and surgical instruments..

- ▶ First-class cleaning power thanks to three enzymes
- ▶ Excellent material compatibility thanks to a neutral pH value
- ▶ Low foaming thanks to non-ionic surfactants
- ▶ Optimal combination of high-quality enzymes (protease = breaks down proteins, lipase = dissolves fatty soils, amylase = removes polysaccharides)
- ▶ Excellent cleaning of all types of thermostable and thermolabile medical devices with excellent material compatibility
- ▶ 5 l canister

307305



DISINFECTANT SOLUTIONS GIGASEPT® FF

ORDER
ONLINE

www.eickemeyer.com

gigasept® FF (new) is a succinic dialdehyde-based preparation for the manual disinfection of surgical and sensitive medical instruments.

- ▶ Free of formaldehyde and GDA
- ▶ Comprehensive microbiological spectrum of activity
- ▶ Excellent material compatibility
- ▶ 7-day shelf life
- ▶ Disinfection of thermolabile and thermostable instruments, e.g., flexible and rigid endoscopes, waterproof ultrasound probes, anesthesia accessories, breathing masks, instruments, etc.

- ▶ 2 l bottle

307315

- ▶ 5 l canister

307320



307315



307320

SCHÜLKE THERMOSEPT® X TRA ENZYME CLEANER

Mildly alkaline high-performance machine cleaner with enzymes.

- ▶ Outstanding cleaning performance thanks to the optimal combination of enzymes and surfactants.
- ▶ All-in-one solution for machine cleaning with a wide range of applications.
- ▶ Designed for thermal reprocessing in washer-disinfectors (WDs) and container washing systems.
- ▶ Highly cost-effective due to low application concentration.
- ▶ Removes organic contaminants such as blood, proteins, tissue residues, mucus, and greasy soils even with low dosages.
- ▶ Offers outstanding material compatibility, even with sensitive materials such as anodized aluminum and non-ferrous metals.
- ▶ Low-foam formulation.
- ▶ Label-free and therefore particularly user-friendly.
- ▶ 5 L canister

307340



SCHÜLKE MIKROZID® UNIVERSAL WIPES GREEN LINE DISINFECTANT WIPES

Low-alcohol rapid disinfection for non-invasive medical devices and surfaces.

- ▶ 100 % plastic-free wipes
- ▶ Multiply sustainable and climate-neutral
- ▶ Maximum material compatibility
- ▶ Effective in seconds
- ▶ Dermatologically tested
- ▶ Bactericidal: 15 sec. – 1 min.
- ▶ Yeastocidal: 1 min.
- ▶ Limited virucidal – 15 sec.

307350

