




Low Lift Exam/Prep Table

veterinary technology for life
Eickemeyer®

Goin' Down



Tie Down Cleat



Low Lift Exam /Prep Table

The hydraulic lift mechanism helps protect your back. The heavy duty frame ensures stability throughout its range. All 4 castors can be locked simultaneously with the foot switch.

- ▶ Minimum height: 49,5 cm
- ▶ Maximum height: 98 cm
- ▶ Maximum load: 100 kg
- ▶ Four lockable castors
- ▶ Table Top dimensions (in cm): L 122 x W 60
- ▶ Includes: 6 tie-down cleats

600200

EICKEMEYER®
Eltastraße 8
78532 Tuttlingen
Germany

T +49 7461 96580 0
F +49 7461 96580 90
E export@eickemeyer.com

TELEPHONE +49 7461 96580 0
www.eickemeyer.com