




Dental Unit
Profident & Profident PLUS

veterinary technology for life
Eickemeyer®

Keep on movin'



174001

Dental Unit Profident & Profident PLUS

The Profident is the result of many years of experience. It offers a wide range of applications for both, general practitioners and dental specialists. It is characterised by its modern design combined with high mobility. The stand is adjustable in height and all handpieces are arranged within easy reach of the operator. When a handpiece is removed from its holder, the Profident automatically supplies air and water to it.

Profident

- ▶ High Speed Turbine (300,000 rpm)
- ▶ Low Speed Air Motor (20,000 rpm) with Straight Handpiece
- ▶ 3 in 1 Syringe (precise finger tip control of water, mist or dry air)
- ▶ Foot Switch (for speed control and water activation)
- ▶ Silent Air Compressor (maintenance free)
- ▶ Includes: Micro Oil Filter

174000

Profident PLUS

- ▶ It offers all the features of the normal Profident (174000). In addition, it includes a built-in EickSonic dental scaler with LED illumination as well as an LED illuminated high speed turbine.

174001

EICKEMEYER®
Eltastraße 8
78532 Tuttlingen
Germany

T +49 7461 96 580 0
F +49 7461 96 580 90
export@eickemeyer.com

TELEPHONE +49 7461 96 580 0
www.eickemeyer.com



Dental Unit Profident & Profident PLUS

veterinary technology for life
Eickemeyer®



174022



174004

Profident Accessories

Low Speed Dental Handpiece with integrated Cooling System

Consisting of:

- ▶ Low Speed Air Motor – With internal water cooling (up to 20,000 rpm, for appropriately configured dental units)
- ▶ Low Speed Straight Handpiece S-Max – 1: 1, with internal water cooling (for appropriately configured dental unit)
- ▶ Low Speed Contra-angle – 1:1, for low speed motor with water cooling (for appropriately configured dental unit)

174022

Starter Kit Profident / Profident PLUS

- ▶ Handpiece Lubricant Oil
- ▶ High Speed Air Motor / Turbine Adapter
- ▶ Low Speed E-Type Adapter
- ▶ 100x Disposable Propfy Heads (HP)
- ▶ Propfy Paste (200 g)
- ▶ Burr Cleaning Brush
- ▶ 2x Dental Burr (HP) Round – Size 23
- ▶ 5x Dental Burr (FG) Round – Size 23
- ▶ 5x Dental Burr (FG) Round Tapered – Size 10
- ▶ 5x Dental Burr (FG) Round Tapered – Size 12
- ▶ 5x Dental Burr (FG) Round Tapered – Size 16

174004

EICKEMEYER®
Eltastraße 8
78532 Tuttlingen
Germany

T +49 7461 96 580 0
F +49 7461 96 580 90
export@eickemeyer.com

TELEPHONE +49 7461 96 580 0
www.eickemeyer.com