



EickView 4-in-1
Endoscopy Monitor

veterinary technology for life
Eickemeyer®

The All-Rounder

EickView 4-in-1 Endoscopy Monitor

The EickView 4-in-1 monitor is a real all-rounder for endoscopy. The device combines a large 22" full HD display, a high-resolution CCD camera, a powerful LED cold light source and digital image data storage in an easy-to-use and robust system.

- ▶ Large 22" full HD screen with 1,920 x 1,080 pixels
- ▶ High-resolution full HD camera with 1,920 x 1,080 pixels
- ▶ Zoom adapter for adjusting the focal length and focus
- ▶ Powerful LED cold light source with 490,000 lux light intensity
- ▶ Integrated image data storage on SD card
- ▶ Recording formats: MP4, JPG
- ▶ HDMI output
- ▶ Incl. remote control
- ▶ Length of fibre optic cable: 2.8 m
- ▶ Dimensions (in cm): W 56 x H 38 x D 15
- ▶ Weight: 15 kg

306105

Camera Cover

- ▶ For single use
- ▶ Telescopic folding
- ▶ Individually sterile packed
- ▶ Dimensions (in cm): L 150 x W 8
- ▶ 50 per pack

306102



306105



306102

EICKEMEYER®
Eltastraße 8
78532 Tuttlingen
Germany

T +49 7461 96580 0
F +49 7461 96580 90
export@eickemeyer.com

TELEPHONE +49 7461 96580 0
www.eickemeyer.com