


Disposable Ring Retractor
with Elastic Hooks

veterinary technology for life
Eickemeyer®

Precision in Surgery



Disposable Ring Retractor with Elastic Hooks

Our ring retractor is effective in its simplicity, ensuring a clear and unobstructed view of the surgical field.

- ▶ **Adaptive Design:** With four flexible hinges, our retractor conforms seamlessly to the patient's body, providing optimal access and comfort during intricate surgeries.
- ▶ **Generous Working Area:** Offers a substantial 14.4 × 14.4 cm field, enabling you to work with greater freedom and precision.
- ▶ **Enhanced Efficiency:** Comes equipped with 8 sharp, elastic hooks, meticulously designed to secure the surgical site swiftly and effectively.
- ▶ **Sterility Assured:** Each retractor is individually packaged sterile, guaranteeing the highest standards of hygiene and patient safety.

162150