



InfusoVet SP Syringe Pump

Fluid Dynamics



InfusoVet SP Syringe Pump

The EICKEMEYER® InfusoVet SP is a highly accurate and reliable syringe pump. The clearly arranged, colour coded buttons make it particularly easy to use.

- ▶ 2.5" backlit LCD display
- ▶ Compatible with standard 5 ml, 10 ml, 20 ml, 30 ml, 50 ml and 60 ml syringes
- ▶ Safety when injecting: precise control of flow rate and injection volume
- ▶ High patient safety through acoustic alarms
- ▶ Adjustable alarm volume
- ▶ Quick start function
- ▶ With clamp for IV pole

Flow rates

- ▶ 5 ml syringe: 0.1 – 150 ml/h
- ▶ 10 ml syringe: 0.1 – 300 ml/h
- ▶ 20 ml syringe: 0.1 – 600 ml/h
- ▶ 30 ml syringe: 0.1 – 900 ml/h
- ▶ 50 ml / 60 ml syringe: 0.1 – 1,500 ml/h

Technical specifications

- ▶ Mains and battery operation, battery life up to 6 hours
- ▶ Splash proof (IP34)
- ▶ Dimensions (in cm): L 33 x D 12.6 x H 11
- ▶ Weight: approx. 1.8 kg

402095