



EickSorber
Anaesthetic Gas Adsorber

veterinary technology for life

Eickemeyer®

For your safety during anaesthesia



213488

EickSorber Anaesthetic Gas Adsorber

- ▶ Activated carbon filter to absorb organic anaesthetic gases
- ▶ Connects to pressure relief valve, thus preventing unwanted gases from entering room atmosphere.
- ▶ Replace filters after exceeding 1,000 g

Single Unit

- ▶ Without exhaust hose

21348801

Single Unit

- ▶ Includes 2 m exhaust hose (Ø 22 mm)

213488

Pack of 6 pieces

- ▶ Includes 2 m exhaust hose (Ø 22 mm)

21348806

EICKEMEYER®
Eltastraße 8
78532 Tuttlingen
Germany

T +49 7461 96 580 0
F +49 7461 96 580 90
export@eickemeyer.com

TELEPHONE +49 7461 96 580 0
www.eickemeyer.com