



HEINE® Ophthalmic
Examination Lamp, 3.5 V

veterinary technology for life
Eickemeyer®

*Look into
my eyes ...*



**HEINE® Ophthalmic
Examination Lamp, 3.5 V**

- ▶ For examination of the anterior segment and cornea in focal light
- ▶ 5 different apertures that help detect erosion, oedema and foreign bodies on the cornea
- ▶ Compatible with HEINE® handles
- ▶ Apertures focused at 100 mm working distance

301180

EICKEMEYER®
Eltastrasse 8
78532 Tuttlingen
Germany

T +49 7461 96580 0
F +49 7461 96580 91
E export@eickemeyer.com

TELEPHONE +49 7461 96580 0
www.eickemeyer.com