




Pawsitioner Positioning Aid

Secure your Position



605925

Pawsitioner Positioning Aid

The Pawsitioner is designed to secure patients for surgery and hands-free diagnostic imaging. It eliminates the need for tie-downs and makes it easy to clip and scrub patients in the prep room then transport them to surgery and back without having to re-position. Each cradle is shaped and contoured to allow for maximum surgical exposure and manoeuvrability while on the surgery table. The Pawsitioner is also ideal for hands-free radiographic positioning, including V-D, left or right lateral, V-D pelvis for OFA evaluation, rostral-caudal skull views and much more without re-positioning the patient for each radiograph. The cradle may also be used to hold patients for laparoscopy, CT, MRI, myelography and ultrasonography.

- ▶ Great insulating characteristics, better exposure and more manoeuvrability
- ▶ Easy to clean

▶ Size XS

- ▶ Length: 60 cm
- ▶ Width: 17 cm

605921

▶ Size M

- ▶ Length: 95 cm
- ▶ Width: 28 cm

605923

▶ Size XL

- ▶ Length: 150 cm
- ▶ Width: 40 cm

605925

▶ Size S

- ▶ Length: 75 cm
- ▶ Width: 22 cm

605922

▶ Size L

- ▶ Length: 115 cm
- ▶ Width: 35 cm

605924

Set of five sizes

605920