



LifeVet CP Capnograph  
and Pulse Oximeter

veterinary technology for life  
**Eickemeyer®**

# Smart Technology made easy



## LifeVet CP Capnograph and Pulse Oximeter

The LifeVet CP handheld Capnograph and Pulse Oximeter measures and displays the expiratory CO<sub>2</sub> concentration as well as the respiratory rate, pulse rate and oxygen saturation in dogs, cats and horses. This compact and portable device is very easy to use with its touchscreen display and intuitive interface.

### Features

- ▶ High resolution 4.3" Touch Screen Display
- ▶ Side-stream Capnography
- ▶ EtCO<sub>2</sub> concentration: 0 – 150 mmHg
- ▶ Respiratory Rate displayed in rpm
- ▶ SpO<sub>2</sub> range: 1 – 100 %
- ▶ Pulse Rate: 25 – 250 bpm
- ▶ Audible and visual alarms with adjustable volume
- ▶ Long-life Lithium battery: Up to 18 h operating time
- ▶ Dimensions (in cm): W 7,7 x H 14,1 x D 3,6
- ▶ Weight: 420 g
- ▶ Includes: CO<sub>2</sub> tube, CO<sub>2</sub> filter, net adapter, USB cable, clip sensor
- ▶ Please note: Please order the exhaust air connection set (321881)

**321880**



### EICKEMEYER®

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**