





EICKEMEYER® I.V. Giving Set

veterinary technology for life
Eickemeyer®

Every drop counts ...



050165



050166

EICKEMEYER® I.V. Giving Set

- ▶ 20 drops/ml
- ▶ Length: 1.50 m
- ▶ Y-connector and filter
- ▶ Includes Luer-Lock connector
- ▶ Pack of 50
- ▶ Individually wrapped sterile

050165

EICKEMEYER® I.V. Giving Set Extension (Heidelberger extension)

- ▶ Length: 1.50 m
- ▶ Individually wrapped sterile
- ▶ Pack of 50

050166

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
E info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk