



DENTO-DORM®
Tracheal Sponge Set

veterinary technology for life
Eickemeyer®

Safe dental treatments



178260

DENTO-DORM® Tracheal Sponge Set

Made from cellular rubber, this is a sealing sponge for the trachea, especially during calculus removal and dental treatments. The DENTO-DORM® Tracheal Sponge has a tightly sealed outlet opening for the tracheal tube. The set includes 48 tracheal sponges and a carrying case.

178260

DENTO-DORM® Tracheal Sponge Refill Pack

- Includes: 12 Coarse Sponges, 3 different sizes

178261

- Includes: 36 Fine Sponges, 6 different sizes

178262



178261



178262

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
E info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk