



ICU for small patients

veterinary technology for life  
**Eickemeyer®**

# *A warm nest*

## *ICU for small patients*

The Vetario T30 unit is an advanced, compact Intensive Care Unit capable of delivering small animal patient warming and recovery. It features full digital temperature control and early warning alarms. The Vetario T30 is ideal for young, sick, injured or weak animals in emergency situations or animal rescue. Its small profile makes it usable in a variety of veterinary settings, from pre-examination warming in consulting rooms, to use in patient homes or mobile ambulances.



- ▶ Digital temperature control, programmable to 0.1 °C
- ▶ Digital temperature reading – choose °C or °F
- ▶ High and low, internal and external temperature alarms
- ▶ Clear, high visibility locking-door
- ▶ Adjustable ventilation
- ▶ Made of anti-microbial plastics
- ▶ Easy to dismantle for thorough cleaning
- ▶ Built-in carry handle
- ▶ 12 V power supply for optional in-car adaptor
- ▶ Stackable with no assembly required
- ▶ Floor area: 240 x 240 mm
- ▶ Total volume: 9.5 litres
- ▶ Dimensions (in mm): H 297 x W 270 x D 275

**219036**

### **EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
F 020 8891 2686  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**