



Digital Cat Scale

veterinary technology for life
Eickemeyer®

Nice and easy



Digital Cat Scale

- ▶ Contoured weighing surface
- ▶ Easy to clean
- ▶ Hold function
- ▶ Tare weighing function
- ▶ LCD display
- ▶ Button-on technology
- ▶ 5 g graduation (0.011 lb / 0.018 oz)
- ▶ 20 kg capacity (44 lb / 704 oz)
- ▶ Automatic switch-off function
- ▶ Toggle between kg / lb / oz
- ▶ Batteries included (2 x 1.5 V)
- ▶ Dimensions (in cm): L 55 x W 31 x H 5.3

909056

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk