




Dr. Mach Examination
Light LED 115

veterinary technology for life
Eickemeyer®

Ray of Light



611150



611151

611152

Dr. Mach Examination Light LED 115

Dr. Mach from Germany has been setting the standards in medical illumination technology for decades. The new generation of lights with LED technology continues this tradition and supports veterinary professionals with superior technology and design.

- ▶ This examination light is ideal for check-ups and minor procedures
- ▶ Flexible gooseneck arm is easy to manoeuvre and allows for accurate light beam positioning
- ▶ Compact design allows its use in small rooms
- ▶ LED technology has a minimum life-span of 40,000 hours
- ▶ Authentic colour rendering to reduce eyestrain
- ▶ Negligible heat production

Additional information:

- ▶ 60,000 lux at 0.5 m distance
- ▶ Light field size: 14 cm at 0.5 m distance
- ▶ Colour Rendering Index R(a) 95
- ▶ Light head: 22 cm
- ▶ Number of LEDs: 7
- ▶ Total power consumption: 10 watts

Stand Model

611150

Wall Model

611151

Table Model

611152

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
E info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk