



Low Lift Scissor Table II

More comfort for your back!



604200

Low Lift Scissor Table II

The EICKEMEYER® Low Lift Scissor Table II is perfect for lifting and transporting larger patients without unnecessary handling and movement. The table is available in 8 different colours and the optional add-ons make this table a real aa-rounder for your practice.

- ▶ Minimal height of 28 cm allowing large dogs to walk onto the table
- ▶ Maximum height: 105 cm
- ▶ Electric operating table helps protect your back
- ▶ Table height adjustable via foot pedal
- ▶ Stainless steel construction
- ▶ Smooth flat top with anti-drip edges and drain hole
- ▶ Wooden plate under table top to absorb vibration and reduce noise
- ▶ Maximum load: 100 kg
- ▶ Tabletop dimensions (in cm): L 130 x W 60
- ▶ Includes: Power Supply and Tie Down Cleats

Variations

Colour	RAL value	Colour name	Item
		Stainless steel	604200
	3000	Flame Red	604201
	6018	Lime Green	604202
	5022	Night Blue	604203
	2001	Red Orange	604204
	6005	Moss Green	604205
	4010	Pink	604206
	9010	Pure White	604207

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk



Low Lift Scissor Table II



604303



604305



604301



604304

Additional Options

In addition to the stationary table, you can choose one or more of the following table options:

Mobile Kit

- ▶ Includes Wheels, Battery and Charger
- ▶ Table not included

604306

Battery Kit

- ▶ For wireless operation for stationary table
- ▶ Sold as battery only

604307

Curved Dental Tray

- ▶ Ergonomic, semi-circular removable drip tray with grid
- ▶ Made of stainless steel
- ▶ With removable tray with grid
- ▶ Tray size (in cm): W 50 x D 35

604303

Integrated Scale Option

- ▶ Weighing range: 0 – 120 kg
- ▶ ATTENTION: Cannot be combined with Tilt Mechanism (item no. 604301)

604305

Electric Tilt Option

- ▶ Foot pedal operated
- ▶ Tilt mechanism up to 35°
- ▶ Cannot be combined with Scale Option (item no. 604305)
- ▶ ATTENTION: Table not included

604301

DIN Table Rail System

- ▶ Lateral
- ▶ Attachable

604304

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk