



EICKEMEYER® Microscope

veterinary technology for life
Eickemeyer®

It's microscope time



EICKEMEYER® Microscope

- ▶ A basic microscope with good optical quality at an affordable price
- ▶ Binocular head Siedentopf type, 30°
- ▶ Widefield eyepieces WF 10x/18 mm
- ▶ Achromatic objectives 4x/0.1, 10x/0.25, 40x/0.65, 100x/1.25 oil
- ▶ Köhlersche illumination system, 12 V/20 W halogen lamp
- ▶ Blue filter
- ▶ Abbe condenser 1.25 N.A.
- ▶ Mechanical stage 190 x 140 mm

710061

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk