



EickFil Ophthalmic Suture Material

veterinary technology for life
Eickemeyer®

Keep Calm and Sew On!



050889

EickFil Ophthalmic Suture Material

EickFil PGA

- Type: Braided and coated multifilament
- Coating: Polycaprolactone and calcium stearate (1 %)
- Tissue reaction: Minimal
- Absorption: The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days. Approximately 50 % of tensile strength remain after 21 days.
- Presentation: Box with 12 sutures

Characteristics:

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure
- Coated multifilament
- Hermetically sealed packing

EickFil PGA

Item No.	USP, Metric, Needle	Colour	Length
050889	- USP 9/0, EP 0.3 3/8 circle, DSP 6.4 mm, 200 µm - Side cutting spatula	violet	30 cm
050890	- USP 6/0, EP 0.7 3/8 circle, DS 12 mm - Reverse cutting	violet	45 cm
050891	- USP 3/0, EP 2 3/8 circle, DS 19 mm easy pass - Reverse cutting	undyed	50 cm
050892	- USP 3/0, EP 2 3/8 circle, DS 24 mm - Reverse cutting	undyed	75 cm
050893	- USP 6/0, EP 0.7 1/2 circle, HR 17 mm - Round bodied taper point	violet	75 cm

Item No.	USP, Metric, Needle	Colour	Length
050894	- USP 10/0, EP 0.2 3/8 circle, 2 x DSP 6.2 mm, 150 µm - Side cutting spatula double armed	violet	30 cm
050895	- USP 9/0, EP 0.3 3/8 circle, 2 x DSP 6.4 mm, 200 µm - Side cutting spatula double armed	violet	30 cm
050896	- USP 8/0, EP 0.4 3/8 circle, 2 x DSP 6 mm, 200 µm - Side cutting spatula double armed	violet	45 cm
050897	- USP 7/0, EP 0.5 3/8 circle, 2 x DSP 8 mm, 200 µm - Side cutting spatula double armed	violet	45 cm
050898	- USP 6/0, EP 0.7 1/4 circle, 2 x VSP 8 mm, 300 µm - Side cutting spatula double armed	violet	45 cm

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk



EickFil Ophthalmic Suture Material

veterinary technology for life
Eickemeyer®



050899



050901

EickFil Nylon

- ▶ Type: Monofilament
- ▶ Composition: Extrusion of polyamide 6.0 or 6.6
- ▶ Coating: None
- ▶ Tissue reaction: Minimal
- ▶ Absorption: Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
- ▶ Presentation: Box with 12 sutures

Characteristics:

- ▶ High tensile strength, non absorbable
- ▶ Synthetic
- ▶ Pliability and softness
- ▶ Monofilament with microscopic uniform diameter
- ▶ Hermetically sealed packing

EickFil Nylon

Item No.	USP, Metric, Needle	Colour	Length
050899	- USP 10/0, EP 0.2 3/8 circle, 2 x DSP 6.4 mm, 220 µm - Side cutting spatula double armed	black	45 cm
050900	- USP 9/0, EP 0.3 3/8 circle, 2 x DSP 6.2 mm, 150 µm - Side cutting spatula double armed	black	45 cm

EickFil Silk

- ▶ Type: Braided multifilament
- ▶ Composition: Braided fibres from the cocoon of the silkworm
- ▶ Coating: Waxes Silicone
- ▶ Tissue reaction: Moderate
- ▶ Absorption: Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.
- ▶ Presentation: Box with 12 sutures

Characteristics:

- ▶ High tensile strength, non absorbable
- ▶ Braided or twisted structure
- ▶ Coated multifilament
- ▶ Hermetically sealed packing

EickFil Silk

Item No.	USP, Metric, Needle	Colour	Length
050901	- USP 4/0, EP 1.5 3/8 circle, DS 19 mm - Reverse cutting	black	45 cm

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk