



Suture: EickFil

veterinary technology for life
Eickemeyer®

Premium suture for the demanding surgeon



050607

Suture: EickFil

- ▶ Braided and Coated Absorbable Suture Material
- ▶ Coating: Calcium Stearate and Polycaprolactone
- ▶ Material: Coated Polyglycolic Acid
- ▶ Absorption by hydrolysis
- ▶ Tensile strength: 84 % at 14 days, 23 % at 28 days
- ▶ Total absorption: 60 to 90 days
- ▶ Colour: Violet
- ▶ 12 packs per box
- ▶ Thread length: 70 cm

Suture: EickFil – Round Bodied			
Article No.	Metric	Length mm	USP
050607	2	26	3/0
050615	2	30	3/0
050608	3	26	2/0
050616	3	30	2/0
050617	3.5	30	0
050618	4	30	1

Suture: EickFil – Reverse Cutting			
Article No.	Metric	Length mm	USP
050650	1	19	5/0
050651	1.5	19	4/0
050658	1.5	24	4/0
050652	2	19	3/0
050659	2	24	3/0
050665	2	30	3/0

Suture: EickFil – Reverse Cutting			
Article No.	Metric	Length mm	USP
050653	3	19	2/0
050660	3	24	2/0
050666	3	30	2/0
050661	3.5	24	0
050667	3.5	30	0
050668	4	30	1

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk