



EickView LED Video
Endoscopes for horses

veterinary technology for life
Eickemeyer®

The Upper Airways



306166



EickView LED Video Endoscopes for horses

The EickView E endoscope range is the perfect solution for examination of the upper respiratory tract in horses. The powerful LED light source and the high-resolution camera of this stand-alone system are integrated in the distal end of the scope. The user-friendly EickView E endoscopes deliver brilliant digital images e.g. of the sinuses or the guttural pouches. The 2.8 mm working channel and the scope length make it ideal for performing biopsies and tracheal washes with ease.

Technical Specifications

- ▶ Diameter: 8.0 mm
- ▶ Working channel: Ø 2.8 mm
- ▶ Depth of field: 5–80 mm
- ▶ Direction/field of view: 0°/110°
- ▶ 3.5" monitor with 3 h battery life
- ▶ SD card for video and image storage

EickView 100E

- ▶ For foals and ponies
- ▶ Working length: 100 cm
- ▶ Bend up/down: 150°/130°

306166

EickView 150E

- ▶ Working length: 150 cm
- ▶ Bend up/down: 130°

306164

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk