



Anaesthesia Masks
for Cats and Dogs

veterinary technology for life
Eickemeyer®

Every breath you take



Anaesthesia Masks for Cats and Dogs

- ▶ Robust, clear perspex
- ▶ Sizes 3–6 supplied with adjustable harness

Complete Set of 6 different sizes

- ▶ Size 1 (Item no. 215401)
OD: 44 mm, ID: 22 mm, depth: 24 mm
- ▶ Size 2 (Item no. 215402)
OD: 57 mm, ID: 30 mm, depth: 29 mm
- ▶ Size 3 (Item no. 215403)
OD: 60 mm, ID: 33 mm, depth: 35 mm
- ▶ Size 4 (Item no. 215404)
OD: 87 mm, ID: 30 mm, depth: 73 mm
- ▶ Size 5 (Item no. 215405)
OD: 110 mm, ID: 42 mm, depth: 97 mm
- ▶ Size 6 (Item no. 215406)
OD: 130 mm, ID: 54 mm, depth: 121 mm

215400

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk