



EICKEMEYER® standard
centrifuge

veterinary technology for life
Eickemeyer®

Spinning Around



EICKEMEYER® Standard Centrifuge

- ▶ Safe and easy to operate
- ▶ Excellent value for money
- ▶ Holds 6 x 15 ml tubes
- ▶ Variable speed: 0 – 4500 RPM
- ▶ Built in timer, running time up to 60 min
- ▶ Safety latch stops the motor when the lid is open

718200

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
E info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk