

ULTRASOUND MAGIC 500 PLUS

A new standard

*High
mobility*

*Compact
shape*



Ultrasound

A NEW STANDARD MAGIC 500 PLUS

We observe the work of the veterinarians in practice and design devices that adapt to the needs of our customers. The MAGIC 500 PLUS meets the highest demands for pregnancy testing in the field.

**High
mobility**

**Compact
shape**

Long Battery Life

Due to the long battery life of the battery supplied, unforeseen extended scanning sessions are not a problem.

Simple Archiving

You can save pictures and videos either directly onto the device or externally via a WLAN connection. The ability to perform measurements on stored images ensures an efficient examination process.

Effective

At only 1.1 kg and an integrated rechargeable battery, it ensures maximum comfort and user friendliness.

“The MAGIC 500 PLUS is ideal for my mobile large animal practice. During the examination in the stable, I have a very good view with the device on my arm.”



ULTRASOUND MAGIC 500 PLUS

Areas of Application

- ▶ Pregnancy testing in cattle and horses

Adjustment Aids for Optimal Work

- ▶ Standard measurements and vet pregnancy calculation
- ▶ Measurements below 1 mm possible (fat layer thickness measurement)
- ▶ Back fat and loin thickness calculation

Special Features

- ▶ Anti-reflective ultrasound image individually adjustable
- ▶ Protected by a splash-proof cover and membrane keypad
- ▶ Wireless connection to an external monitor
- ▶ WLAN connection to PC/laptop/tablet/smartphone
- ▶ Incl. 6.5 MHz Linear Rectal Probe



Technical Specifications

- | | | |
|---|--|---|
| ▶ 5.9" LED monitor | ▶ Higher diagnostic reliability with Cine Memory (256 frames) | ▶ Dimensions (in cm): W 16 x H 16 x D 7 |
| ▶ Images and video clips can be saved on the device | ▶ Incl. glove for left hand (right glove is an optional extra) | ▶ Weight: 1.1 kg |
| ▶ Unlimited storage of images via WLAN connection | ▶ Unit can be easily hung around your neck via mounting rings and straps | 303790 |

Ultrasound

ACCESSORIES FOR MAGIC 500 PLUS

*Transducer technology with outstanding image quality.
Our probes are available in several designs so we can
offer you a solution that best meets your requirements.*



PROBE MICRO CONVEX

► 6.5 MHz


303791



PROBE LINEAR T-SHAPED

► 7.5 MHz

303793



PROBE LINEAR RECTAL

► Included

► Frequency range: 5.5 – 9.0 MHz

303792



GLOVE (LEFT HAND) FOR MAGIC 500 PLUS


► Included

30379002



GLOVE (RIGHT HAND) FOR MAGIC 500 PLUS

30379001



Examination of a pregnant mare