



HEINE GAMMA®
Stethoscopes

veterinary technology for life
Eickemeyer®

*I can feel your
heart beat...*

MADE IN
GERMANY
5 YEARS
WARRANTY



310054

HEINE GAMMA® Stethoscopes

For excellent acoustic performance.

- ▶ Combined double chest piece with two different frequency ranges
- ▶ High-quality membrane for excellent acoustic performance and sound transmission
- ▶ Conical acoustic head on the other side
- ▶ Single-tube system
- ▶ Latex-free

HEINE GAMMA® Stethoscope, large

- ▶ Membrane size: Ø 44 mm
- ▶ Cone size: Ø 30 mm

310054

HEINE GAMMA® Stethoscope, small

- ▶ Membrane size: Ø 31 mm
- ▶ Cone size: Ø 30 mm

310055



310055



310056

HEINE GAMMA® Cardio Stethoscope

Exceptional acoustic performance in the cardiac frequency range.

- ▶ Combined double chest piece with double membrane on one side and a conical acoustic head on the other
- ▶ High-performance double membrane, 45 mm diameter for perfect skin contact for exceptional sound reproduction
- ▶ Two-in-one tube system for enhanced performance
- ▶ Latex-free

310056



EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
E info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk